



# Chapter 10

## PLASmid TAXonomic PCR (PlasTax-PCR), a Multiplex Relaxase MOB Typing to Assort Plasmids into Taxonomic Units

Raquel Cuartas, Teresa M. Coque, Fernando de la Cruz, and M. Pilar Garcillán-Barcia

### Abstract


Plasmids transmissible by conjugation are responsible for disseminating antibiotic-resistance genes, making plasmid detection relevant for pathogen tracking. We describe the use of a multiplex PCR method for the experimental identification of specific plasmid taxonomic units (PTUs) of transmissible plasmids. The PCR primers were designed to target conserved segments of the relaxase MOB gene of PTUs encoding adaptive traits for enterobacteria (antimicrobial resistance, virulence, and metabolism). In this way, PlasTax-PCR detects the presence of these plasmids and allows their direct assignment to a PTU.

**Key words** Plasmid taxonomic units, Bacterial conjugation, Horizontal gene transfer, Relaxase, MOB family, Plasmid typing

---

### 1 Introduction

Plasmids are critical vehicles in disseminating antimicrobial resistance (AMR) [1]. Thus, their detection and classification are crucial for molecular epidemiology to track AMR beyond the boundaries of specific bacterial clones. Plasmids transmissible by conjugation have a distinctive characteristic, they encode a MOB relaxase that recognizes a cognate plasmid sequence, the *nic* site within the origin of transfer (*oriT*) [2]. Nine MOB relaxase classes are currently described [3]. Still, just five of them comprise more than 95% of conjugative relaxases present in plasmids hosted in the order Enterobacterales [4]. To detect transmissible plasmids, a PCR-based method, Degenerate-Primer MOB typing (DPMT), was developed [5, 6] (Fig. 1). Applied to collections of clinical and environmental isolates of enterobacteria [8–12], DPMT detected plasmids with relaxase sequences identical or nonidentical



**Fig. 1** DPMT and PlasTax-PCR scheme concepts. Phylogenetic tree representations of the five relaxase MOB classes prevalent in Enterobacterales are schematized. Families into each MOB class with members in Enterobacterales are represented by colored triangles. Those identified by the DPMT Scheme [6, 7] are indicated in the left panel. At the right panel, the same families are depicted, but only the tips contain colored triangles corresponding to the PTUs included in the PlasTax-PCR scheme

to those previously known and classified them into broad MOB families.

Plasmids have been recently assorted into taxonomic units (PTUs), which gather members with a common genomic backbone

[4]. One thousand seven hundred and seventy out of the 2535 plasmids hosted in the order Enterobacterales were included in 83 PTUs. Out of them, 55 PTUs did not correspond to any known incompatibility (Inc) group. Nevertheless, 50 out of these 83 PTUs were MOB<sup>+</sup>, according to MOBscan [3]. Besides, each PTU was characterized by a single relaxase MOB type, while replication functions within a given PTU showed considerable variation [4]. These facts endorse the use of the MOB relaxase sequences to classify transmissible plasmids into their corresponding PTUs. Here, we describe a PCR-based method, Plasmid Taxonomic PCR (PlasTax-PCR), to detect plasmid relaxases from PTUs circulating in clinical enterobacteria. This method targets specifically the MOB relaxases corresponding to 19 PTUs, and of two additional groups without PTU assignment (IncT and ColA-like colicin plasmids) (Fig. 1). For each group, the coding sequences of the N-terminal relaxase domain (5' 900 nucleotides) were aligned and conserved specific blocks were chosen to design the primers. As proof of principle, this method was applied to detect transmissible plasmids in a series of *E. coli* ST131 clinical isolates, a polyclonal cluster able to harbor a high variability of plasmids (Fig. 2).

---


## 2 Materials

### 2.1 Solutions (See Notes 1 and 2)

1. InstaGene™ Matrix (Bio-Rad).
2. DNA polymerase and its corresponding reaction buffer, e.g., BioTaq™ DNA polymerase (Bioline) (see Note 3).
3. 50 mM MgCl<sub>2</sub>.
4. 10 mM dNTP.
5. PCR-grade H<sub>2</sub>O (see Note 4).
6. 100 μM solutions for each primer. PlasTax-PCR primers are listed in Table 1 (see Notes 5 and 6).
7. Agarose basic, DNase free.
8. 10× TBE (Tris/Borate/EDTA) buffer, pH 8.3.
9. Intercalating agent, e.g., SafeView™ Classic (ABM) (5 μL stock/100 mL gel).
10. Ladder for DNA electrophoresis, e.g., HyperLadder™ 1 kb (Bioline).
11. DNA loading buffer: 30% (v/v) glycerol, 0.25% (w/v) bromophenol blue.

### 2.2 Equipment

12. UV-visible spectrophotometer, e.g., NanoDrop™ 2000 (Thermo Scientific).
13. A thermocycler to carry out the PCR reactions.



**Fig. 2** Analysis of PlasTax-PCR reactions. For each multiplex, a 2% agarose gel is shown. M: HyperLadder™ 1Kb. Positive controls (see Table 1): individual plasmids, indicated above the corresponding lane, and a mix of the individual plasmids in the same PCR sample (mix C+). Negative control: no DNA. Samples: PCRs from clinical *E. coli* ST131 isolates, whose genomes are fully sequenced (strains FV9873, E35BA, E2022, E61BA, and FV13998 from [14]) or not (*E. coli* isolates from [15, 16])

- 14. An agarose electrophoresis system.
- 15. A gel imaging system to visualize fluorescent nucleic acid stains.

### 3 Methods

#### 3.1 Preparation of the PlasTax-PCR DNA Template

Total DNA was extracted from *Escherichia coli* cultures bearing the plasmids listed in the second column of Table 1 using InstaGene™ Matrix (Bio-Rad) and following the manufacturer's recommendations (see Note 7).

1. Centrifuge 50 µL of a saturated bacterial culture at 13,800 rcf for 1 min.

Table 1  
Targets, primers, and conditions for PlasTax-PCR

Multiplex	MOB family/ PTU/plasmid prototype	Forward primer (5' → 3')	Reverse primer (5' → 3')	Amplicon size (bp)	T <sub>ann</sub> (°C)
1	MOB <sub>P11</sub> /PTU-P1/R751	CTGGCGAACCAGCACACGA	GTAGTCGCGCTCCAGGGCC	297	62
	MOB <sub>P13</sub> /PTU-L/M/pCTX-M3	AAATCATTTAGGGGCCGAAGGG	CATACTGGTGATCGGACATGCC	436	
	MOB <sub>P12</sub> /PTU-II/R64	CATGAAGGACGGCTGCGAATG	ATTTCCGTCATCTGAGACAGG	629	
	MOB <sub>P3</sub> /PTU-X1/pOLA52	AGTCATCAGAAAATGGTCGTAAG TCAGCT	CTTTCCCGCGAACATAAAGTCCC TACA	823	
	16S rRNA <sup>a</sup>	CTCCTACGGGAGGCAGCAG	GACGGCGGTGTGTRCA	1034	
2	MOB <sub>F11</sub> /PTU-W/R388	CATAGGGCGGGCTGCAAGC	GTGGCCCTCGCCGATATTCC	169	62
	MOB <sub>H12</sub> /PTU-C/pIP1202	GTATGCCGCGTTTGTTCACCTG	GGCCACCGGATAAAAGTAACG	340	
	MOB <sub>F12</sub> /PTU-F <sub>E</sub> /R100 <sup>b</sup>	GTTGAGTTTCAGCGTTGTAAA TCGG	AGCACGTTCTCGGAAAAATCCG	567	
	MOB <sub>H11</sub> /PTU-HI2/R478	ATCTCAGGAGAATGATGCAACCTC TG	AAGACATACCGGGTTTAGGA TTCCG	712	
	16S rRNA <sup>a</sup>	CTCCTACGGGAGGCAGCAG	GACGGCGGTGTGTRCA	1034	
3	MOB <sub>P12</sub> /PTU-B/O/K/Z/R387	CGGACAAAGCTGTTTTTCCCGTA	TCCTGGCTCCAGGAAAGCCAG	253	62
	MOB <sub>H11</sub> /PTU-HI1A/R27	TGAGCAGTCTATTCTTTTCGTTG TACAGAGAC	GGCGCCCATGTTAATGTTTCACTC	405	
	MOB <sub>F12</sub> /PTU-F <sub>K1</sub> - PTU-F <sub>K2</sub> / pKpQIL <sup>c</sup>	GGATAAGGCGCTGTTTACGGAAC TG	GGTCTGCAGCCCTCCCTGACC	600	
	16S rRNA <sup>a</sup>	CTCCTACGGGAGGCAGCAG	GACGGCGGTGTGTRCA	1034	
	MOB <sub>P51</sub> /-/ColA; pCK02	GGGGCGAAACATCCGACT	ACCGTTCTGTCTCGCTGAAT	142	
4	MOB <sub>C11</sub> /PTU-E4/ColDF13	CGCACCGGCTGGCCG	TGACCTGCGTCGCCCGG	270	64
	MOB <sub>Q41</sub> /PTU-E10/pIGWZ12	CCCCTGCCCTGGTGTACGAACC	TATGAACGATGGTCTATCTCTCC TGATAAC	394	
	16S rRNA <sup>a</sup>	CTCCTACGGGAGGCAGCAG	GACGGCGGTGTGTRCA	1034	

(continued)

Table 1  
(continued)

Multiplex	MOB family/PTU/plasmid prototype	Forward primer (5' → 3')	Reverse primer (5' → 3')	Amplicon size (bp)	T <sub>ann</sub> (°C)
5	MOB <sub>Q4</sub> /PTU-E20/pIGMS5	CCCGTGAGATTTCGCCGCGAG	CTTTCCTGCGAATCCCGTTTCC	176	65
	MOB <sub>F11</sub> /PTU-N1/R46	CTTGATATAACCACGATTACCCGCC	CTTCACGCACAGACGGC	355	
	MOB <sub>H12</sub> /-/Rts1	AGCTGAATTGGGTTTTCGGC	CCGCTCAATCTTTTCAGTTTCGG	460	
	16S rRNA <sup>a</sup>	CTCCTACGGGAGGCAGCAG	GACGGCGGTGTGTRCA	1034	
6	MOB <sub>p51</sub> /PTU-E1/ColE1	AAAAGGCCGTCAGGATGTGATTCA	GCGTCTCTGATTTTCGTCTCG TTTG	280	64
	MOB <sub>F12</sub> /PTU-F <sub>8</sub> /pSLT	GGCGGCAACAAACACCCGC	AACCTTCTCTTTCAGCACCGCG	411	
	MOB <sub>F12</sub> /PTU-F <sub>9</sub> /pMT	CGGACTCAGGACGGGGCG	CCACTCGGCCATGCGCTG	630	
	16S rRNA <sup>a</sup>	CTCCTACGGGAGGCAGCAG	GACGGCGGTGTGTRCA	1034	
7	MOB <sub>F12</sub> /PTU-F <sub>6</sub> /F <sup>b</sup>	GCTGGGCAGCATGGGAGAAC	CAATCTGATTAGCGTACACATTC TCAATG	509	62
	16S rRNA <sup>a</sup>	CTCCTACGGGAGGCAGCAG	GACGGCGGTGTGTRCA	1034	

<sup>a</sup>16S rDNA primers are 357F and 1391R, described by [13]

<sup>b</sup>PTU-F<sub>6</sub> presents low conservation in the gene repertoire of its members [4], and thus two different primer pairs were designed to target this group

<sup>c</sup>This is the only case in which a PlasTax-PCR primer pair targets plasmids from two different PTUs. Plasmids from the related group PTU-F<sub>K3</sub> are not targeted

2. Add a volume of 200  $\mu\text{L}$  of InstaGene™ Matrix to the pellet.
3. Incubate at 56 °C for 15–30 min.
4. Vortex at high speed for 10 s.
5. Incubate at 100 °C for 8 min.
6. Repeat **step 4**.
7. Centrifuge at 13,800 rcf for 3 minutes and recover the supernatant.
8. Repeat **step 7** to eliminate traces of the matrix.
9. Quantify the DNA (*see Note 8*).
10. Store the supernatant at  $-20\text{ }^{\circ}\text{C}$  (*see Note 9*).

### **3.2 PlasTax-PCR for Enterobacterial Plasmids**

1. **Preparation of the reaction mixture for the multiplex PCRs.** For a reaction of 50  $\mu\text{L}$ , add 100 ng of genomic DNA (*see Notes 8 and 10*), 5  $\mu\text{L}$  of 10 $\times$  reaction buffer, 1.5  $\mu\text{L}$  of 50 mM  $\text{MgCl}_2$  (final concentration 1.5 mM), 1  $\mu\text{L}$  of 10 mM dNTP (final concentration 0.2 mM), 0.5  $\mu\text{L}$  of 100  $\mu\text{M}$  primers (final concentration 1  $\mu\text{M}$ ), 1 U of BioTaq polymerase (Bio-line), and ddH<sub>2</sub>O up to 50  $\mu\text{L}$  (*see Note 4*).
2. **PCR running** (*see Note 11*). The amplification program includes an initial denaturation step 5 min at 94 °C, followed by 30 cycles of 30 s at 94 °C + 30 s at annealing temperature (*see Table 1*) + 30 s at 72 °C, and a final extension at 72 °C during 7 min.
3. **Visualization of the amplicons.** Load 10  $\mu\text{L}$  of the PCR reaction and 2  $\mu\text{L}$  of the DNA loading buffer onto a 2% agarose gel containing an intercalating agent. Load at least one lane of DNA ladder in each gel (*see Note 12*). Separate by electrophoresis at room temperature at voltage 100 V for 40 min. Visualize under ultraviolet light (*see Note 13*).

---

## **4 Notes**

1. Pipette all reagents and samples using filter tips. The preparation of the reactions before the DNA amplification should be ideally carried out in a clean area (pre-PCR) to reduce the chance of contamination.
2. Aliquot reagents to avoid multiple freeze-thaw cycles and the contamination of master stocks.
3. No need to use extremely low error-rate DNA polymerases if amplicons are not going to be sequenced.
4. The use of commercially available nuclease-free water is preferred.

5. A thorough process of primer design is recommendable. A nucleotide multiple sequence alignment of the relaxase gene portion encoding the N-terminal relaxase domain is the basis to detect high identity stretches for primer designing. It is useful to design several primer pairs to detect a single PTU, which render amplicons of different sizes, in order to ease their combination in multiplex reactions. Sequences that produce intra- or inter-oligonucleotide secondary structures should be avoided. The melting temperature ( $T_m$ ) of the primers should be as close as possible, favoring  $T_m > 60^\circ\text{C}$  to avoid nonspecific targets. The primer pair selected should have good performance in both the individual and the multiplex PCR reaction.
6. The individual amplicons included in a multiplex reaction should differ in size, ideally 150–200 bp, as to be clearly distinguished in a 2% agarose gel. Take it into account when designing the primers.
7. A purification method rendering PCR-quality template DNA is advised, instead of picking and adding a bacterial colony directly into the PCR master mix.
8. We prepared aliquots of template DNA at 20 ng/ $\mu\text{L}$  to add the same volume to each PCR reaction.
9. DNAs extracted using InstaGene Matrix can be stored at  $-20^\circ\text{C}$  for a month without degradation.
10. When a mix of template DNAs is used in the reaction, 100 ng of each DNA sample is added.
11. Positive controls for each reaction should be used in parallel. A list of GenBank accession numbers of plasmids that can be targeted by each primer pair and can thus be used as positive controls, according to the criteria explained in the footnote to Table 2, is provided:

*For multiplex 1:*

PTU-P1: NC\_013176.1, NC\_005088.1, NC\_017908.2,  
 NC\_019263.1, NC\_019264.1, NC\_019283.1,  
 NC\_016968.1, NC\_016978.1, NC\_004956.1,  
 NC\_006830.1, NC\_021077.1, NC\_019312.1,  
 NC\_001735.4, NC\_024998.1, NC\_007353.2,  
 NC\_010935.1, NC\_019320.1, NZ\_CP017760.1,  
 NC\_020994.1, NZ\_CP015373.1, NZ\_CP009797.1,  
 NZ\_CP014846.1, NZ\_CP019238.1, NC\_014911.1,  
 NC\_008766.1, NC\_014641.1, NC\_007337.1,  
 NZ\_CP021650.1, NC\_008385.1.

PTU-L/M: NC\_004464.2, NC\_005246.1,  
 NC\_011641.1, NC\_019063.1, NC\_019154.1,  
 NC\_019344.1, NC\_019346.1, NC\_019368.1,  
 NC\_019889.1, NC\_021078.1, NC\_021488.1,  
 NC\_023027.1, NC\_024997.1, NC\_025134.1,



**Table 2**  
**PlasTax-PCR efficiency for each targeted PTU**

Multiplex	PTU	Number of PTU members (RefSeq84) <sup>a</sup>	Number of targeted PTU members <sup>b</sup>
1	P1	32	29
	L/M	58	56
	I1	120	119
	X1	49	24
2	W	4	4
	C	92	87
	F <sub>E</sub>	198	52
	HI2	33	29
3	B/O/K/Z	20	14
	HI1A	8	8
	F <sub>K1</sub> - F <sub>K2</sub>	68/32	37/29
4	E4	40	39
	E10	23	21
5	E20	9	8
	N1	60	60
6	E1	51	45
	F <sub>S</sub>	56	43
	F <sub>Y</sub>	38	10
7	F <sub>E</sub>	198	94

<sup>a</sup>Data are taken from Supplementary Table S4 [4]  
<sup>b</sup>Number of plasmids that are potentially targeted by the corresponding primers. Take into account that not all plasmids from a PTU contain a MOB relaxase. Plasmids are considered targeted if the primers potentially anneal to their sequence producing a fragment with the size indicated in Table 1. Hybridization mismatches are allowed, except if they occur in the 3' 12 nucleotides of the primers, in which case the plasmid was ruled out. According to these criteria, a list of plasmids that can be used as positive controls is provided in **Note 11**

NZ\_CP007733.1, NZ\_CP009852.1, NZ\_CP009857.1,  
NZ\_CP010365.1, NZ\_CP011593.1, NZ\_CP011599.1,  
NZ\_CP011609.1, NZ\_CP011614.1, NZ\_CP011632.1,  
NZ\_CP011640.1, NZ\_CP014298.1, NZ\_CP014698.2,  
NZ\_CP015071.1, NZ\_CP015075.2, NZ\_CP016927.1,  
NZ\_CP017282.1, NZ\_CP017288.1, NZ\_CP017853.1,  
NZ\_CP017932.1, NZ\_CP017936.1, NZ\_CP018315.1,  
NZ\_CP018342.1, NZ\_CP018449.1, NZ\_CP018452.1,  
NZ\_CP018461.1, NZ\_CP018669.1, NZ\_CP018690.1,  
NZ\_CP018700.1, NZ\_CP018706.1, NZ\_CP018712.1,  
NZ\_CP018717.1, NZ\_CP018723.1, NZ\_CP018736.1,  
NZ\_CP018974.1, NZ\_CP019841.1, NZ\_CP020844.1,  
NZ\_CP021742.1, NZ\_CP022147.1, NZ\_CP022150.1,  
NZ\_CP022153.1, NZ\_CP022826.1, NZ\_KX118608.1.

PTU-II: NC\_013120.1, NC\_014383.1, NC\_019044.1,  
 NC\_019043.1, NC\_019061.1, NC\_019097.1,  
 NC\_002122.1, NC\_022885.1, NC\_023326.1,  
 NC\_024980.1, NC\_023915.1, NC\_025140.1,  
 NC\_025144.1, NC\_025147.1, NC\_025180.1,  
 NC\_025198.1, NC\_025142.1, NC\_025143.1,  
 NC\_025176.1, NC\_024975.1, NC\_024976.1,  
 NC\_024977.1, NC\_024979.1, NC\_024978.1,  
 NC\_024955.2, NC\_019123.1, NC\_019131.1,  
 NC\_019137.1, NC\_019099.1, NC\_005014.1,  
 NC\_015965.1, NC\_019111.1, NC\_023899.1,  
 NC\_023900.1, NC\_019104.1, NC\_023329.1,  
 NC\_022267.1, NC\_023275.1, NC\_023276.1,  
 NC\_023290.1, NC\_020991.1, NC\_022742.1,  
 NC\_018659.1, NZ\_CP019220.1, NC\_011419.1,  
 NC\_016904.1, NC\_017637.1, NC\_017665.1,  
 NC\_017642.1, NZ\_CP009580.1, NZ\_CP006641.1,  
 NZ\_CP010317.1, NZ\_CP012627.1, NZ\_CP015161.1,  
 NZ\_CP015916.1, NZ\_CP015996.1, NZ\_CP015838.1,  
 NZ\_CP018116.1, NZ\_CP018122.1, NZ\_CP018110.1,  
 NZ\_CP010130.1, NZ\_CP010233.1, NZ\_CP019215.1,  
 NC\_017675.1, NC\_011081.1, NC\_017718.1,  
 NC\_011077.1, NC\_021811.1, NC\_021813.2,  
 NZ\_CP009566.1, NZ\_CP012039.1, NZ\_CP012835.1,  
 NZ\_LN890525.1, NZ\_CP012923.1, NZ\_CP012936.1,  
 NZ\_CP013224.1, NZ\_CP012929.1, NZ\_CP013221.1,  
 NZ\_CP014662.1, NZ\_CP016516.1, NZ\_CP016520.1,  
 NZ\_CP016533.1, NZ\_CP016522.1, NZ\_CP016568.1,  
 NZ\_CP016572.1, NZ\_CP016585.1, NZ\_CP016387.1,  
 NZ\_CP016407.1, NZ\_CP016411.1, NZ\_CP016413.1,  
 NZ\_CP016409.1, NZ\_CP019205.1, NZ\_CP019207.1,  
 NZ\_CP010142.1, NZ\_CP019173.1, NZ\_CP018946.1,  
 NZ\_CP018975.1, NZ\_CP018993.1, NZ\_CP014972.2,  
 NZ\_CP014622.1, NZ\_KX443694.1, NZ\_LT838199.1,  
 NC\_032100.1, NZ\_CP018774.2, NZ\_CP021208.1,  
 NZ\_CP021533.1, NZ\_CP021693.1, NZ\_CP021739.1,  
 NZ\_CP021841.1, NZ\_CP021845.1, NZ\_CP021882.1,  
 NZ\_CP018625.1, NZ\_CP020494.1, NZ\_CP019691.1,  
 NZ\_CP022456.1, NZ\_CP014096.1, NZ\_CP010831.1,  
 NZ\_CP014494.1, NZ\_CP007651.1.

PTU-XI: NC\_010378.1, NC\_010421.1, NC\_010422.1,  
 NC\_010860.1, NC\_011204.1, NC\_011739.1,  
 NC\_013503.1, NC\_015472.1, NC\_016036.1,  
 NC\_019013.1, NC\_019046.1, NC\_019067.1,  
 NC\_019088.1, NC\_019096.1, NC\_019106.1,  
 NC\_019256.1, NC\_024961.1, NZ\_CP011431.1,  
 NZ\_CP012734.1, NZ\_CP014974.1, NZ\_CP019180.1,  
 NZ\_CP020088.1, NZ\_CP020341.1, NZ\_CP020836.1.

*For multiplex 2:*

PTU-W: NC\_010643.1, NC\_009982.1, NC\_010716.1,  
NC\_028464.1.

PTU-C: NC\_008612.1, NC\_008613.1, NC\_009139.1,  
NC\_009140.1, NC\_012690.1, NC\_012692.1,  
NC\_012693.1, NC\_016974.1, NC\_016976.1,  
NC\_017645.1, NC\_018994.1, NC\_019045.2,  
NC\_019065.1, NC\_019066.1, NC\_019069.1,  
NC\_019107.1, NC\_019116.1, NC\_019118.1,  
NC\_019121.1, NC\_019153.1, NC\_019158.1,  
NC\_019375.1, NC\_019380.1, NC\_020180.1,  
NC\_021667.1, NC\_021815.1, NC\_022372.1,  
NC\_022377.1, NC\_022522.2, NC\_022652.1,  
NC\_023291.1, NC\_023898.1, NC\_023908.1,  
NZ\_CP003998.1, NZ\_CP006661.1, NZ\_CP007486.1,  
NZ\_CP007636.1, NZ\_CP008790.1, NZ\_CP009409.2,  
NZ\_CP009560.1, NZ\_CP009562.1, NZ\_CP009564.1,  
NZ\_CP009567.1, NZ\_CP009570.1, NZ\_CP009868.1,  
NZ\_CP010373.2, NZ\_CP010391.1, NZ\_CP011429.1,  
NZ\_CP011540.1, NZ\_CP011622.1, NZ\_CP012682.1,  
NZ\_CP013324.1, NZ\_CP014295.1, NZ\_CP014658.1,  
NZ\_CP014775.1, NZ\_CP014978.1, NZ\_CP015139.1,  
NZ\_CP015394.1, NZ\_CP015835.1, NZ\_CP016013.1,  
NZ\_CP016036.1, NZ\_CP017055.1, NZ\_CP017987.1,  
NZ\_CP018318.1, NZ\_CP018689.1, NZ\_CP018698.1,  
NZ\_CP018704.1, NZ\_CP018710.1, NZ\_CP018716.1,  
NZ\_CP018722.1, NZ\_CP018817.1, NZ\_CP018956.1,  
NZ\_CP019441.1, NZ\_CP020049.1, NZ\_CP020056.1,  
NZ\_CP021206.1, NZ\_CP021551.1, NZ\_CP021709.1,  
NZ\_CP021719.1, NZ\_CP021835.1, NZ\_CP021853.1,  
NZ\_CP021936.1, NZ\_CP021952.1, NZ\_CP021956.1,  
NZ\_CP022126.1, NZ\_CP022359.1, NZ\_LT904892.1.

PTU-Fe: NC\_019424.1, NC\_019095.1, NC\_019090.1,  
NC\_019072.1, NC\_019071.1, NC\_019057.1,  
NC\_018998.1, NC\_017630.1, NC\_016039.1,  
NC\_013727.1, NC\_013655.1, NC\_013542.1,  
NC\_013175.1, NC\_011812.1, NC\_011749.1,  
NC\_009133.1, NC\_007941.1, NC\_005327.1,  
NC\_002134.1, NZ\_CP018125.1, NZ\_CP018119.1,  
NZ\_CP018113.1, NZ\_CP018107.1, NZ\_CP014496.1,  
NZ\_CP017287.1, NZ\_CP015072.1, NZ\_CP014523.1,  
NC\_008460.1, NZ\_CP020340.1, NZ\_CP014493.1,  
NZ\_CP012139.1, NZ\_CP021938.1, NZ\_CP021871.1,  
NZ\_CP021180.1, NZ\_CP020117.1, NZ\_CP019009.1,  
NZ\_CP018982.1, NZ\_CP018954.1, NZ\_CP018952.1,  
NZ\_CP016035.1, NZ\_CP015140.1, NZ\_CP015077.1,  
NZ\_CP015070.1, NZ\_CP010882.1, NZ\_CP009860.1,  
NZ\_CP009579.1, NZ\_CP008715.1, NC\_013951.1,

NC\_025177.1, NZ\_CP021690.1, NZ\_CP014498.1,  
NZ\_CP010239.1.

PTU-HI2: NC\_005211.1, NC\_009838.1,  
NC\_010870.1, NC\_012555.1, NC\_012556.1,  
NC\_019114.1, NC\_024983.1, NZ\_CP008825.1,  
NZ\_CP008899.1, NZ\_CP008906.1, NZ\_CP011062.1,  
NZ\_CP011601.1, NZ\_CP012170.1, NZ\_CP012931.1,  
NZ\_CP015833.1, NZ\_CP016526.1, NZ\_CP016764.1,  
NZ\_CP016838.1, NZ\_CP019214.1, NZ\_CP019443.1,  
NZ\_CP019559.1, NZ\_CP019647.1, NZ\_CP020493.1,  
NZ\_CP021177.1, NZ\_CP021209.1, NZ\_CP022165.1,  
NZ\_CP022533.1, NZ\_CP022696.1, NZ\_CP023143.1.

*For multiplex 3:*

PTU-B/O/K/Z: NZ\_CP023144.1, NZ\_CP015141.1,  
NZ\_CP013024.1, NZ\_CP009107.1, NC\_025138.1,  
NC\_022996.1, NC\_022992.1, NC\_022371.1,  
NC\_018995.1, NC\_014843.1, NC\_007365.1,  
NZ\_CP005999.1, NC\_011754.1, NZ\_CP018772.1.

PTU-HIIA: NC\_002305.1, NC\_003384.1,  
NC\_009981.1, NC\_013365.1, NC\_016825.1,  
NC\_023289.2, NZ\_CP022495.1, NZ\_LT904879.1.

PTU-F<sub>K1</sub> - F<sub>K2</sub>: NC\_019390.1, NC\_024992.1,  
NC\_020132.1, NC\_021654.1, NC\_009649.1,  
NC\_022078.1, NZ\_CP007729.1, NZ\_CP008800.1,  
NZ\_CP008829.1, NZ\_CP008930.1, NZ\_CP009777.1,  
NZ\_CP010393.1, NZ\_CP010574.1, NZ\_CP011577.1,  
NZ\_CP011977.1, NZ\_CP011990.1, NZ\_CP015386.1,  
NZ\_CP015823.1, NZ\_CP018355.1, NZ\_CP018365.1,  
NZ\_CP018424.1, NZ\_CP018430.1, NZ\_CP018434.1,  
NZ\_CP018441.1, NZ\_CP018460.1, NZ\_CP018693.1,  
NZ\_CP019773.1, NZ\_CP020072.1, NZ\_CP020109.1,  
NZ\_CP020838.1, NZ\_CP020842.1, NZ\_CP021540.1,  
NZ\_CP021544.1, NZ\_CP021713.1, NZ\_CP021752.1,  
NZ\_CP021834.1, NZ\_CP022145.1, NC\_009650.1,  
NZ\_CP015395.1, NZ\_LT216438.1, NZ\_CP022698.1,  
NZ\_CP022693.1, NZ\_CP019774.1, NZ\_CP018992.1,  
NZ\_CP018989.1, NZ\_CP015824.1, NZ\_CP014765.1,  
NZ\_CP014669.1, NZ\_CP014650.1, NZ\_CP011991.1,  
NZ\_CP011986.1, NZ\_CP009875.1, NZ\_CP009115.1,  
NZ\_CP008833.1, NZ\_CP008830.1, NZ\_CP007730.1,  
NC\_025187.1, NC\_025167.1, NC\_025166.1,  
NC\_023906.1, NC\_023905.1, NC\_023904.1,  
NC\_023903.1, NC\_019165.1, NC\_019155.1,  
NC\_014016.1.

*For multiplex 4:*

MOB<sub>P51</sub>/-/ColA, pCK02: NC\_009794.1,  
NC\_001373.1, NC\_016151.1.

PTU-E4: NC\_002119.1, NC\_009793.1, NC\_018953.1,  
 NC\_019156.1, NC\_019159.1, NC\_020182.1,  
 NC\_021666.1, NC\_022083.1, NZ\_CP006928.1,  
 NZ\_CP007728.1, NZ\_CP008828.1, NZ\_CP008832.1,  
 NZ\_CP009772.1, NZ\_CP009778.1, NZ\_CP009873.1,  
 NZ\_CP009877.1, NZ\_CP010394.1, NZ\_CP011979.1,  
 NZ\_CP011982.1, NZ\_CP011984.1, NZ\_CP011988.1,  
 NZ\_CP011993.1, NZ\_CP014648.1, NZ\_CP015384.1,  
 NZ\_CP015391.1, NZ\_CP018353.1, NZ\_CP018425.1,  
 NZ\_CP018431.1, NZ\_CP018433.1, NZ\_CP018445.1,  
 NZ\_CP019776.1, NZ\_CP020112.1, NZ\_CP020840.1,  
 NZ\_CP021543.1, NZ\_CP021717.1, NZ\_CP021779.1,  
 NZ\_CP021837.1, NZ\_CP022576.1, NZ\_LT216440.1.

PTU-E10: NZ\_CP016041.1, NZ\_CP019024.1,  
 NZ\_CP018996.1, NZ\_CP018972.1, NZ\_CP018942.1,  
 NC\_019134.1, NZ\_CP011432.1, NZ\_CP011140.1,  
 NZ\_CP018980.1, NZ\_CP019648.1, NZ\_CP015142.1,  
 NZ\_CP012629.1, NC\_010885.1, NZ\_CP018966.1,  
 NZ\_CP018208.1, NZ\_CP010878.1, NZ\_CP006634.1,  
 NC\_011411.1, NZ\_HG941720.1, NC\_010486.1,  
 NZ\_CP018941.1.

*For multiplex 5:*

PTU-E20: NC\_010883.1, NC\_011977.1,  
 NC\_012882.1, NC\_020412.1, NC\_022585.1,  
 NC\_023325.1, NZ\_CP014098.1, NC\_013367.1.

PTU-N1: NC\_003292.1, NC\_007682.3, NC\_009131.1,  
 NC\_009132.1, NC\_009980.1, NC\_011383.1,  
 NC\_011385.1, NC\_011617.1, NC\_014208.1,  
 NC\_014231.1, NC\_014368.1, NC\_015599.1,  
 NC\_019033.1, NC\_019082.1, NC\_019087.1,  
 NC\_019098.1, NC\_019124.1, NC\_019888.1,  
 NC\_020086.1, NC\_020088.1, NC\_021622.1,  
 NC\_021660.2, NC\_021664.2, NC\_022374.1,  
 NC\_022375.1, NC\_023909.1, NC\_023910.1,  
 NC\_024967.1, NC\_024974.1, NC\_025019.1,  
 NC\_025183.1, NC\_025186.1, NC\_032101.1,  
 NZ\_CP008901.1, NZ\_CP008908.1, NZ\_CP009853.1,  
 NZ\_CP009858.1, NZ\_CP009862.1, NZ\_CP009864.1,  
 NZ\_CP009867.1, NZ\_CP009874.1, NZ\_CP009881.1,  
 NZ\_CP011589.1, NZ\_CP014524.1, NZ\_CP017725.1,  
 NZ\_CP018442.1, NZ\_CP018945.1, NZ\_CP018959.1,  
 NZ\_CP018963.1, NZ\_CP018977.1, NZ\_CP019006.1,  
 NZ\_CP019026.1, NZ\_CP020059.1, NZ\_CP020119.1,  
 NZ\_CP021899.1, NZ\_KX062091.1, NZ\_KX154765.1,  
 NZ\_KX276209.1, NZ\_KX397572.1, NZ\_LT838197.1.

MOB<sub>H12</sub>/-/Rts1: NC\_003905.1.

*For multiplex 6:*

PTU-E1: NC\_001371.1, NC\_002809.1, NC\_003079.1,  
 NC\_003457.1, NC\_005019.1, NC\_005970.1,  
 NC\_008488.1, NC\_009791.1, NC\_010485.1,  
 NC\_010672.1, NC\_011214.1, NC\_011407.1,  
 NC\_011799.1, NC\_013363.1, NC\_014235.1,  
 NC\_017321.1, NC\_017636.1, NC\_017654.1,  
 NC\_017655.1, NC\_017661.1, NC\_017662.1,  
 NC\_017721.1, NC\_018997.1, NC\_019076.1,  
 NC\_019078.1, NC\_019102.1, NC\_019136.1,  
 NC\_019250.1, NC\_019357.1, NC\_019982.1,  
 NC\_020251.1, NC\_025004.1, NC\_025026.1,  
 NC\_025178.1, NZ\_CP010175.1, NZ\_CP012628.1,  
 NZ\_CP014198.2, NZ\_CP015917.1, NZ\_CP016039.1,  
 NZ\_CP016513.1, NZ\_CP016519.1, NZ\_CP016584.1,  
 NZ\_CP017845.1, NZ\_CP018943.1, NZ\_CP019025.1.

PTU-F<sub>S</sub>: NC\_002638.1, NC\_003277.2, NC\_006855.1,  
 NC\_007208.1, NC\_012124.1, NC\_013437.1,  
 NC\_014476.2, NC\_016855.1, NC\_016858.1,  
 NC\_016861.1, NC\_016864.1, NC\_017054.1,  
 NC\_017720.1, NC\_019001.1, NC\_019108.1,  
 NC\_019109.1, NC\_019112.1, NC\_022570.1,  
 NZ\_AP014566.1, NZ\_CP007489.1, NZ\_CP007582.1,  
 NZ\_CP008745.1, NZ\_CP012345.1, NZ\_CP012348.1,  
 NZ\_CP013721.1, NZ\_CP014050.1, NZ\_CP014357.1,  
 NZ\_CP014359.1, NZ\_CP014537.1, NZ\_CP014577.1,  
 NZ\_CP014968.1, NZ\_CP014970.1, NZ\_CP014973.1,  
 NZ\_CP014976.1, NZ\_CP014980.1, NZ\_CP015158.1,  
 NZ\_CP016390.1, NZ\_CP017618.1, NZ\_CP017729.1,  
 NZ\_CP020923.1, NZ\_CP022137.1, NZ\_LN999012.1,  
 NZ\_LT855377.1.

PTU-F<sub>Y</sub>: NC\_006323.1, NC\_009378.1,  
 NZ\_CP006746.1, NZ\_CP006749.1, NZ\_CP006752.1,  
 NZ\_CP006756.1, NZ\_CP006760.1, NZ\_CP006779.1,  
 NZ\_CP009714.1, NZ\_CP010248.1.

*For multiplex 7:*

PTU-F<sub>E</sub>: NZ\_LT838198.1, NZ\_CP017726.1,  
 NZ\_CP021847.1, NZ\_CP021843.1, NZ\_CP018994.1,  
 NZ\_CP014489.1, NZ\_CP012636.1, NC\_025106.1,  
 NC\_019089.1, NZ\_CP010223.1, NZ\_CP010215.1,  
 NZ\_CP019007.1, NZ\_CP018455.1, NZ\_CP017981.1,  
 NZ\_CP017633.1, NZ\_CP014273.1, NZ\_CP014271.1,  
 NZ\_CP011496.1, NZ\_CP010158.1, NC\_025175.1,  
 NC\_025139.1, NC\_024956.1, NC\_019073.1,  
 NC\_004998.1, NC\_002483.1, NZ\_CP012113.1,  
 NZ\_CP010316.1, NZ\_CP010192.1, NC\_013362.1,  
 NC\_011747.1, NZ\_CP013026.1, NC\_010409.1,  
 NZ\_CP015914.1, NZ\_CP012683.1, NZ\_CP012626.1,  
 NZ\_CP012928.1, NZ\_CP011916.1, NZ\_CP006001.1,

NC\_018966.1, NC\_018954.1, NC\_017627.1,  
 NC\_014382.1, NC\_011076.1, NC\_006671.1,  
 NZ\_CP021880.1, NZ\_CP021733.1, NZ\_CP021289.1,  
 NZ\_CP021204.1, NZ\_CP019028.1, NZ\_CP018990.1,  
 NZ\_CP018958.1, NZ\_CP016498.1, NZ\_CP010372.1,  
 NZ\_CP010232.1, NC\_023315.1, NC\_020271.1,  
 NC\_014615.1, NC\_010558.1, NC\_010488.1,  
 NZ\_CP011019.1, NZ\_CP010184.1, NZ\_CP010181.1,  
 NZ\_CP019561.1, NZ\_CP021737.1, NZ\_CP011064.1,  
 NZ\_CP018978.1, NZ\_CP013833.1, NZ\_CP013832.1,  
 NZ\_CP011493.1, NZ\_CP010138.1, NC\_019037.1,  
 NC\_014384.1, NC\_007675.1, NZ\_CP012632.1,  
 NZ\_CP013027.1, NC\_019122.1, NC\_019117.1,  
 NZ\_CP022166.1, NZ\_CP020934.1, NZ\_CP019909.1,  
 NZ\_CP019018.1, NZ\_CP018775.2, NZ\_CP015239.1,  
 NZ\_CP010141.1, NZ\_CP010123.1, NZ\_CP006635.1,  
 NZ\_CP005931.1, NC\_025179.1, NC\_017659.1,  
 NC\_012944.1, NC\_009837.1, NC\_010720.1,  
 NZ\_CP009167.1, NC\_017640.1.

12. For gels with more than ten lanes, include the DNA ladder in several lanes at different positions.
13. The amplicons obtained in a multiplex PCR can be confirmed by repeating the corresponding simplex PCR. Its product can be purified and sequenced by standard Sanger sequencing to verify the result.

---

## Acknowledgments

This work was financed by the Spanish Ministry of Science and Innovation (PID2020-117923GB-I00) to FdIC and MPG-B, and by the European Commission, the Joint Programming Initiative in Antimicrobial Resistance and the Ministry of Science and Innovation (UE/STARCS AC16/00039, UE/ST131TS AC16/00043, and PI18/01942) to TMC.

## References

1. de la Cruz F, Davies J (2000) Horizontal gene transfer and the origin of species: lessons from bacteria. *Trends Microbiol* 8(3):128–133. <http://www.ncbi.nlm.nih.gov/pubmed/10707066>
2. Garcillán-Barcia MP, Francia MV, de la Cruz F (2009) The diversity of conjugative relaxases and its application in plasmid classification. *FEMS Microbiol Rev* 33(3):657–687. <http://www.ncbi.nlm.nih.gov/pubmed/19396961>
3. Garcillán-Barcia MP, Redondo-Salvo S, Vielva L, de la Cruz F (2020) MOBscan: automated annotation of MOB relaxases. In: de la Cruz F (ed) *Horizontal gene transfer: methods and protocols*, Methods in molecular biology. Science+Business Media, LLC, part of Springer Nature, New York, NY, pp 295–308
4. Redondo-Salvo S, Fernández-López R, Ruiz R, Vielva L, de Toro M, Rocha EPC et al (2020) Pathways for horizontal gene transfer in bacteria revealed by a global map of their plasmids. *Nat Commun* 11(1):3602. <http://www.ncbi.nlm.nih.gov/pubmed/32681114>

5. Garcillán-Barcia MP, Alvarado A, de la Cruz F (2011) Identification of bacterial plasmids based on mobility and plasmid population biology. *FEMS Microbiol Rev* 35(5):936–956. <http://www.ncbi.nlm.nih.gov/pubmed/21711366>
6. Alvarado A, Garcillán-Barcia MP, de la Cruz F (2012) A degenerate primer MOB typing (DPMT) method to classify gamma-proteobacterial plasmids in clinical and environmental settings. *PLoS One* 7(7):e40438. <http://www.ncbi.nlm.nih.gov/pubmed/22792321>
7. Villa L, Carattoli A (2020) Plasmid typing and classification. Springer, New York, pp 309–321. [https://doi.org/10.1007/978-1-4939-9877-7\\_22](https://doi.org/10.1007/978-1-4939-9877-7_22)
8. Valverde A, Cantón R, Garcillán-Barcia MP, Novais A, Galán JC, Alvarado A et al (2009) Spread of bla(CTX-M-14) is driven mainly by IncK plasmids disseminated among *Escherichia coli* phylogroups A, B1, and D in Spain. *Antimicrob Agents Chemother* 53(12):5204–5212. <http://www.ncbi.nlm.nih.gov/pubmed/19786598>
9. Mata C, Miró E, Alvarado A, Garcillán-Barcia MP, Toleman M, Walsh TR et al (2012) Plasmid typing and genetic context of AmpC  $\beta$ -lactamases in Enterobacteriaceae lacking inducible chromosomal ampC genes: findings from a Spanish hospital 1999–2007. *J Antimicrob Chemother* 67(1):115–122. <http://www.ncbi.nlm.nih.gov/pubmed/21980067>
10. Coelho A, Piedra-Carrasco N, Bartolomé R, Quintero-Zarate JN, Larrosa N, Cornejo-Sánchez T et al (2012) Role of IncHI2 plasmids harbouring blaVIM-1, blaCTX-M-9, aac(6')-Ib and qnrA genes in the spread of multiresistant *Enterobacter cloacae* and *Klebsiella pneumoniae* strains in different units at Hospital Vall d'Hebron, Barcelona, Spain. *Int J Antimicrob Agents* 39(6):514–517. <http://www.ncbi.nlm.nih.gov/pubmed/22481058>
11. Garcillán-Barcia MP, Ruiz del Castillo B, Alvarado A, de la Cruz F, Martínez-Martínez L (2015) Degenerate primer MOB typing of multiresistant clinical isolates of *E. coli* uncovers new plasmid backbones. *Plasmid* 77:17–27. <http://www.ncbi.nlm.nih.gov/pu>
12. Adelowo OO, Caucci S, Banjo OA, Nnanna OC, Awotipe EO, Peters FB et al (2018) Extended Spectrum Beta-lactamase (ESBL)-producing bacteria isolated from hospital wastewaters, rivers and aquaculture sources in Nigeria. *Environ Sci Pollut Res* 25(3):2744–2755. <https://doi.org/10.1007/s11356-017-0686-7>
13. Turner S, Pryer KM, Miao VPW, Palmer JD (1999) Investigating deep phylogenetic relationships among cyanobacteria and plastids by small subunit rRNA sequence analysis. *J Eukaryot Microbiol* 46(4):327–338. <https://doi.org/10.1111/j.1550-7408.1999.tb04612.x>
14. Lanza VF, de Toro M, Garcillán-Barcia MP, Mora A, Blanco J, Coque TM et al (2014) Plasmid flux in *Escherichia coli* ST131 sublineages, analyzed by plasmid constellation network (PLACNET), a new method for plasmid reconstruction from whole genome sequences. *PLoS Genet* 10(12):e1004766
15. Coque TM, Novais A, Carattoli A, Poirel L, Pitout J, Peixe L et al (2008) Dissemination of clonally related *Escherichia coli* strains expressing extended-Spectrum  $\beta$ -lactamase CTX-M-15. *Emerg Infect Dis* 14(2):195–200. [http://wwwnc.cdc.gov/eid/article/14/2/07-0350\\_article.htm](http://wwwnc.cdc.gov/eid/article/14/2/07-0350_article.htm)
16. Novais A, Pires J, Ferreira H, Costa L, Montenegro C, Vuotto C et al (2012) Characterization of globally spread *Escherichia coli* ST131 isolates (1991 to 2010). *Antimicrob Agents Chemother* 56(7):3973–3976. <http://www.ncbi.nlm.nih.gov/pubmed/22491693>